

Certificate of Conformity

Indicator Solution

Universal Indicator for pH Determination of Aqueous Solutions (pH4 - 10)

Product No:	UN101	Date of Measurement:	13/01/2023
Lot No:	UN10123A1	Date of Sample Receipt:	13/01/2023
Expiry Date*:	28/01/2025	Date of Manufacture*:	13/01/2023

Sensitivity Test:

Colour Change: pH 4 Pink / pH 5 Orange / pH 6 Yellow / pH 7 Green / pH 8 Turquoise / pH 9 Blue / pH 10 Purple.

Compliant

Products are manufactured under NSAI registered I.S EN ISO9001:2015 Quality Systems, registration no: 19.2769

Date of Issue of the Certificate :
17/01/2023

Senior QC Technician
MALONE Eileen

Eileen Malone

All raw materials used to prepare this product are of high purity.

*The detail above is based on information supplied in writing by Reagecon Manufacturing.

Tested by Reagecon Quality Control Laboratories for Reagecon Manufacturing

This Certificate must not be reproduced except in full. Rev-QL001