

Certificate of Analysis

Analytical Volumetric Solutions

Sodium Oxalate 0.05M (0.1N)

Product No:	NAX0051
Lot No:	NAX0520E1
Expiry Date:	28/05/2022
Mean Molarity:	0.04988 M @ 20°C
Date of Measurement:	27/05/2020

Specification:
0.04975 - 0.05025M

Method:
Standardised in accordance with in-house potentiometric titration.

Traceability:
This volumetric solution was checked by means of Potassium Permanganate Analytical Volumetric Standard. This volumetric standard is directly traceable to a Standard Reference Material of National Institute of Standards and Technology (USA), 136f Potassium Dichromate. All raw materials used to prepare this product are of high purity.

Reference:
Standard Methods of Preparation, Standardisation and Storage of Standard Solutions for Chemical Analysis ASTM E200-08.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 27/05/2020

Eileen Malone

QC Technician Eileen Malone

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-20E14

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com