

Certificate of Analysis

Redox Standard

124 ± 5mV @25°C (Pt v Ag-AgCl Ref 3M KCl)

Product No:	RS124
Lot No:	RS12421D1
Expiry date:	28/04/2023

Redox potential:	129 mV @ 25°C
Date of measurement:	30/04/2021

Measurement Limits:
119 - 129mV @25°C

Measurement Cell:
Pt/Sample/3MKCl/AgCl/Ag

Method:
Direct potentiometry
Temperature versus redox value

T°C	mV
10	154
15	143
20	136
25	124
30	115

Reference: Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 30/04/2021

A handwritten signature in blue ink, appearing to read 'Orla Sheehan'.

QC Technician Orla Sheehan

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21D13