

CERTIFIKÁT O ANALÝZE

Produkt : **KYSELINA CHLORISTÁ ANALPURE®**
 Objednací číslo : SAc11002
 Číslo šarže : 7019
 Exspirace : 04/04/2023 (dd/mm/yyyy)

 Barevnost : < 10 APHA
 Koncentrace : 68 %

Stopové nečistoty (ppb):

Analyt	Maximální koncentrace	Aktuální hodnota	Analyt	Maximální koncentrace	Aktuální hodnota
Aluminium	1	< 0.5	Neodymium	0.5	< 0.5
Antimony	0.5	< 0.1	Nickel	1	< 0.1
Arsenic	0.5	< 0.5	Niobium	Pouze informativně	< 0.5
Barium	1	< 0.1	Palladium	0.5	< 0.5
Beryllium	0.5	< 0.1	Platinum	0.5	< 0.5
Bismuth	0.5	< 0.1	Potassium	1	< 0.5
Cadmium	1	< 0.1	Praseodymium	0.5	< 0.5
Calcium	1	< 1	Rhodium	0.5	< 0.5
Cerium	0.5	< 0.5	Rubidium	0.5	< 0.5
Cesium	0.5	< 0.5	Samarium	0.5	< 0.5
Cobalt	0.5	< 0.1	Scandium	0.5	< 0.5
Copper	0.5	< 0.5	Silver	1	< 0.1
Dysprosium	0.5	< 0.5	Sodium	1	< 0.5
Erbium	0.5	< 0.5	Strontium	0.5	< 0.1
Europium	0.5	< 0.5	Tantalum	Pouze informativně	< 1
Gadolinium	0.5	< 0.5	Tellurium	0.5	< 0.5
Gallium	0.5	< 0.5	Terbium	0.5	< 0.5
Gold	0.5	< 0.5	Thallium	0.5	< 0.5
Hafnium	Pouze informativně	< 0.5	Thorium	1	< 0.1
Holmium	0.5	< 0.5	Thulium	0.5	< 0.5
Indium	0.5	< 0.5	Tin	1	< 0.5
Iron	1	< 0.5	Titanium	1	< 0.1
Lanthanum	0.5	< 0.5	Tungsten	Pouze informativně	< 0.5
Lead	1	< 0.1	Uranium	0.5	< 0.1
Lithium	0.5	< 0.1	Vanadium	0.5	< 0.5
Lutetium	0.5	< 0.5	Ytterbium	0.5	< 0.5
Magnesium	1	< 0.5	Yttrium	0.5	< 0.5
Manganese	1	< 0.1	Zinc	1	< 0.5
Molybdenum	0.5	< 0.1	Zirconium	0.5	< 0.1



Datum: **04. 04. 2021**

Odpovědná osoba: **Vladimír Patera**