



Certificate of Analysis

Chemical Oxygen Demand - COD Reagents

COD Calibration Standard 20ppm

Product No:
Lot No:
Expiry date * :

COD20
COD2021K1
28/10/2023

Mean Assay:
Date of Measurement:
Date of Sample Receipt:
Date of Manufacturing * :

20.10 ppm as COD
19/10/2021
19/10/2021
19/10/2021

* The detail above is based on information supplied in writing by Reagecon Manufacturing.

Specification:

19.80-20.20ppm

Method:

Standardised in accordance to in-house method, acidimetric.

Traceability:

This volumetric solution was checked by means of Sodium Hydroxide Analytical Volumetric Standard. This volumetric standard is directly traceable to a Standard Reference Material of National Institute of Standards and Technology, NIST SRM 84L, Potassium Hydrogen Phthalate.

Reference:

Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Date of Issue of the Certificate:

27/10/2021

Eileen Malone

QC Technician **Eileen Malone**

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21K17

Tested by Reagecon Quality Control
Laboratories for Reagecon Manufacturing

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com