

Certificate of Analysis

Analytical Volumetric Solution

Potassium Hydroxide 0.1M (0.1N) in Ethanol

Product No:
Lot No:
Expiry Date:

ETKOH01F
EK20121G1
28/07/2022

Mean Molarity:
Date of Measurement:
Specification:
0.0996 - 0.1004M

0.1000 M @ 20°C
03/08/2021

Method:

Tested by acidimetric titration, to in-house method TPATB.

Traceability:

This volumetric solution was checked by means of Hydrochloric Acid Analytical Volumetric Standard. This volumetric standard is directly traceable to a Standard Reference Material of National Institute of Standards and Technology (USA), 723e Tris(hydroxymethyl) aminomethane. All raw materials used to prepare this product are of high purity.

Date of issue of the certificate:

03/08/2021

Eileen Malone

QC Technician Eileen Malone

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1. This Certificate must not be reproduced except in full. Rev-21H5

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com