

Certificate of Analysis

Redox Standard

600± 5mV @25°C (Pt v Ag-AgCl Ref 3M KCl)

Product No:
Lot No:
Expiry date:

RS600
RS60021D1
28/04/2023

Redox potential:
Date of measurement:
Measurement Limits:
595 - 605mV @ 25°C

597 mV @ 25°C
08/04/2021

Measurement Cell:
Pt/Sample/3MKCl/AgCl/Ag

Method:
Direct potentiometry
Temperature versus redox value:
T°C mV
10 577
15 585
20 593
25 600
30 609

Reference:Vogels Textbook of Quantitative Inorganic Analysis. Fifth Edition.

Products are manufactured under an NSAI registered I.S EN ISO9001:2015 Quality System, registration no: 19.2769

Date of issue of the certificate: 08/04/2021

A handwritten signature in blue ink, appearing to read 'Reguymen'.

QC Technician Patricia de Guzman

Complementary information relative to this product is available at www.reagecon.com
Page 1 of 1.This Certificate must not be reproduced except in full. Rev-21C65

Reagecon Diagnostics Ltd.
Shannon Free Zone, Shannon, Co. Clare, Ireland.
Tel +353 61 472622, Fax: +353 61 472642
Email: sales@reagecon.ie, www.reagecon.com