



Shannon Free Zone,
Shannon, Co. Clare, Ireland
Tel: +353 61 472622
Fax: +353 61 472642
Email: sales@reagecon.ie
www.reagecon.com

CERTIFICATE OF ANALYSIS

PRODUCT: 25 % Brix Standard (Stabilised)
PRODUCT NO: BS25S
LOT NO: RIBS2521M1
DATE OF TEST: 15th December 2021
EXPIRY DATE*: 28th December 2022
DATE OF MANUFACTURING*: 15th December 2021
DATE OF SAMPLE RECEIPT: 15th December 2021
SPECIFICATION: 24.80 – 25.20% Brix @ 20°C
MEAN: 24.99% Brix @ 20°C

Refractive index @ 20°C is 1.37231 using sodium light at 589.3nm.

TEST METHOD:

This certificate relates solely to the sample as received by the laboratory, bearing the product code and lot number given above, taken at time of manufacture. The Brix value of the standard was determined using a high performance calibrated, temperature controlled refractometer according to TPRIA 01.

This test result is traceable to high purity ISO 17034 Certified Reference Materials.

This certificate relates solely to the sample as received by the laboratory, bearing the product code and lot number given above. The uncertainty of measurement has been calculated not to exceed $\pm 0.11\%$ Brix value (20°C) at 95% confidence level, i.e. coverage factor $k=2$.

ACCREDITATION:

Reagecon Diagnostics Ltd. is accredited by the Irish National Accreditation Board, under scope 264T, for the test method, TPRIA 01, used to generate the above result. This accreditation is intended only to certify that Reagecon has the Quality Management Systems in place to ensure that each individual test result generated using TPRIA 01 is technically valid and is supported by appropriate uncertainty measurements.

INSTRUCTIONS FOR USE:

- Prior to use, the contents of the dropper bottle should be brought to room temperature and mixed thoroughly by inversion to ensure a homogeneous mix.
- Remove cap and discharge 4-5 droplets of liquid to waste to clear the nozzle.
- Replace cap immediately after use.
- Product must be stored at 2-8°C when not in use.
- DO NOT FREEZE

This Product is Stabilised by the addition of 0.01% of preservative.

*The detail above is based on information supplied in writing by Reagecon Manufacturing.

Tested by Reagecon Quality Control Laboratories for Reagecon Manufacturing.

Approved By: QC Technician



Date: 15th December 2021

This certificate must not be reproduced except in full.